


NEW LAB
SOLUTIONS



SDLE - SPME DUAL LAYER EXTRACTION

PRODUCT FEATURES

- Easy installation above the standard heated agitator
- Temperature set by Chronos software within the sample preparation method
- Special vial cap to better heat transfer surface
- Up to - 20 °C (this specific are defined by the thermo exchanger installed)



SRA 
INSTRUMENTS
ANALYTICAL SOLUTIONS

www.srainstruments.it

 **CHROMLINE**

SRA Instruments Spa - IT
Ph.: +39 02.92143.58
info@srainstruments.com

SRA Instruments Sas - FR
Tel: +33 4.78442947
info@sra-instruments.com